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WITH THE GUARANTEE OF
TELIC GROUP

EN. BITE BLOCK FOR ENDOTRACHEAL TUBES AND LARYNGEAL MASKS

INTENDED USE

Accessory to prevent pressure on the endotracheal tube/probe due to biting, in case of oral intubation. The product is commercialized under the following references:

[REF]	Ø (mm)	Tube / Probe	Patient
7600	19	CH/FR 8-9-10	Adult
7650	15	CH/FR 6-7-8	Paediatric

INSTRUCTIONS FOR USE

- Secure the endotracheal tube according to medical standard protocols.
- Pull the pilot balloon outside the mouth, keeping the pilot balloon tube in a straight line together with the endotracheal tube.
- Place the MORDEDOR-MO bite block next to the endotracheal tube and the pilot balloon tube. Be secure that the bite block does not kink or press the pilot balloon tube, especially at the insertion point to the endotracheal tube.
- Use plastic medical adhesive tape to fix the bite block to the endotracheal tube. Excessive pressure when taping may distort the endotracheal tube, especially in case of paediatric tubes.
- Take note of the position of the bite block in relationship with the endotracheal tube, to control that no migration of the bite block is occurring.
- When removing the endotracheal tube, first of all, detach the bite block and dispose according to hospital standard protocols for contaminated waste, in accordance with regulations in force.

USE TIME

Replace the product after 72h of intubation for hygienic purposes.

WARNINGS AND PRECAUTIONS

- This product should only be used by or under the direction of a health care professional familiar with its correct placement and use. Any allergic or other skin reactions should be assessed by a doctor immediately.
- Only use on intact, clean skin (do not use on open wounds, lesions or infected or inflamed areas).
- Not for internal use. In case of accidental ingestion of the product, seek immediate medical attention.
- Check the cuff pressure regularly. If the cuff is not filling correctly, check that the bite block is neither pushing against nor bending the pilot balloon tube.
- If the bite block needs replacing or repositioning, check that the pilot balloon tube is in good condition, especially at its point of insertion in the endotracheal tube.
- Inspect the bite block regularly. Check that the bite block has not become displaced and is neither pushing against nor bending the pilot balloon tube.
- Check that the bite block is not damaging the patient's lips. The bite block should be positioned as close to the centre of the patient's mouth as possible.
- In patients with filled teeth, bridges or dentures, the benefits should be weighed up against the possibility of permanent damage to the teeth. Introduce and remove the bite block carefully.
- Single-use product. Do not reuse.
- Non-sterile product. This product's sterilisation is not guaranteed.
- The reuse or sterilisation of this product could cause:
 - Cross infections due to the presence of biological residues;
 - Changes in the materials; or
 - Loss of product functionality.
- Do not use if the individual packaging is open or damaged.
- After checking package integrity, check the contents. If the packaging is damaged or has visible defects, do not use the product and return it to TELIC, S.A.U.
- Check the expiry date on the packaging. Do not use after the expiry date.

DISPOSAL

The product must be disposed according to hospital protocols for contaminated waste, in accordance with regulations in force. Dispose the package according to recyclable residues management.

WARRANTY

TELIC, S.A.U. guarantees that the product complies with the Regulation (EU) 2017/745 on Medical devices and was manufactured following the procedures of the TELIC, S.A.U. quality management system certified under ISO standard 13485. TELIC, S.A.U. accepts no liability for any medical costs or direct or indirect harm due to the lack of operation or incorrect operation of the products when used without observing the instructions for use. We advise you to duly report any operational problem or product fault to the Quality Assurance Department of TELIC, S.A.U. In the event of any serious incident related to the product, this must be reported immediately to TELIC, S.A.U. (incidents@telic.es) and/or the competent authority of the Member State in which the incident occurs.

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